



Recombinant Human hepatitis A virus genotype IA Genome polyprotein

Product Code	CSB-MP318397HJF
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P13672
Product Type	Recombinant Protein
Immunogen Species	Human hepatitis A virus genotype IA (isolate LCDC-1) (HHAV) (Human hepatitis A virus (isolate Human/China/LCDC-1/1984))
Purity	≥85% (SDS-PAGE)
Sequence	QVGDDSGGFS TTVSTEQNVP DPQVGIKGA NRGKMDVSGV QAPVGAITTI EDPVLAKKVP ETFPELKPGE SRHTSDHMSI YKFMGRSHFL CTFTFNSNNK EYTFPITLSS TSNPPHGLPS TLRWFFNLQ LYRGPLDLTI IITGATDVDG MAWFTPVGLA VDPWVEKAS ALSIDYKTAL GAVRFNTRRT GNIQIRLPWY SYLYAVSGAL DGLGDKTDST FGLVSIQIAN YNHSDEYLSF SCYLSVTEQS EFYFPRAPLN SNAMLSTEM MSRIAAGDLE SSVDDPRSEE DRRFESHIEC RKPKYKELRLE VGKQRLKYAQ EELSNEVLPP PRKMKGLFSQ S
Source	Mammalian cell
Protein Names	Recommended name: Genome polyprotein Cleaved into the following 4 chains: 1. Protein VP3 Alternative name(s): P1C Virion protein 3 Protein VP1-2A Alternative name(s): PX Protein VP1 Alternative name(s): P1D Virion prote
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.