



Recombinant 30S ribosomal protein S8P (rps8p)

Product Code	CSB-EP318463MEZ
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P14038
Product Type	Recombinant Protein
Immunogen Species	Methanococcus vanniellii
Purity	>85% (SDS-PAGE)
Sequence	MSLMDPLANA LNHVSNCEGV GKNVAYLKPA SKLIGRVLKV MQDQGYIGNF EYIEDGKAGV YKVDLIGQIN KCGAVKPRYA VKYQEFEKFE KRYLPAKGF LLIVSTPKGL MTHDEARTAG VGGRLISYVY
Source	E.coli
Target Names	rps8
Protein Names	Recommended name: 30S ribosomal protein S8P
Expression Region	1-130
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.