



Recombinant L-2-hydroxyisocaproate dehydrogenase

Product Code	CSB-BP318502LAO
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P14295
Product Type	Recombinant Protein
Immunogen Species	Weissella confusa (Lactobacillus confusus)
Purity	≥85% (SDS-PAGE)
Sequence	ARKIGIIGL GNVGAAVAHG LIAQGVADDY VFIDANEAKV KADQIDFQDA MANLEAHGNI VINDWAALAD ADVVISTLGN IKLQQDNPTG DRFAELKFTS SMVQSVGTNL KESGFHGV LV VISNPVDVIT ALFQHVTGFP AHKVIGTGTL LDTARMQRAV GEAFDLDP RS VSGYNLGEHG NSQFVAWSTV RVMGQPIVTL ADAGDIDLAA IEEEEARKGGF TVLNGKGYTS YGVATSAIRI AKAVMADAHA ELVVSNNRRDD MGYLSYP AI IGRDGVLAET TLDLTTDEQE KLLQSRDYIQ QRFDEIVDTL
Source	Baculovirus
Protein Names	Recommended name: L-2-hydroxyisocaproate dehydrogenase Short name= L-HicDH EC= 1.1.1.-
Expression Region	2-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.