



Recombinant *Oryza sativa* subsp. *japonica* Glutelin type-B 1 (GluB1-A)

Product Code	CSB-EP318506OFG
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P14323
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	QLFNPS TNPWHSPRQG SFRECRFDRL QAFEPLRKVR SEAGVTEYFD EKNELFQCTG TFVIRRVIQP QGLLVPRYTN IPGVVYIIQG RGSMLTFPG CPATYQQQFQ QFSSQGQSQS QKFRDEHQKI HQFRQGDIVA LPAGVAHWFY NDGDPIVAV YVYDVNNAN QLEPRQKEFL LAGNNRAQQ QQVYGSSIEQ HSGQNIFSGF GVEMLSEALG INAVAAKRLQ SQNDQRGEII HVKNGLQLLK PTLTQQQEQA QAQDQYQQVQ YSERQQTSSR WN
Source	E.coli
Target Names	GluB1-A
Protein Names	Recommended name: Glutelin type-B 1 Cleaved into the following 2 chains: 1. Glutelin type-B 1 acidic chain 2. Glutelin type-B 1 basic chain
Expression Region	25-302
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.