



# Recombinant Human Catalase (CAT)

<b>Product Code</b>	CSB-BP004563HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P04040
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>ADSRDPASD QMQHWKEQRA A QKADVLTTG AGNPVGDKLN VITVGPRGPL  LVQDVVFTDE MAHFDRERIP ERVVHAKGAG AFGYFEVTHD ITKYSKAKVF  EHIGKKTPIA VRFSTVAGES GSADTVRDPR GFAVKFYTED GNWDLVGNNT  PIFFIRDPIF FPSFIHSQKR NPQTHLKDPD MVWDFWVSLRP ESLHQVSFLF  SDRGIPDGHR HMNGYGSHTF KLVNANGEAV YCKFHYKTDQ GIKNLSVEDA  ARLSQEDPDY GIRDLFNAIA TGKYPVSWTFY IQVMTFNQAE TFPFNPFDLT  KVWPHKDYPL IPVGKLVVLR NPVNYFAEVE QIAFDPSNMP PGIEASPKDM  LQGRLFAYPD THRHRLGPNY LHIPVNCYPYR ARVANYQRDG  PMCMQDNQGG APNYYPNVNSFG APEQQPSALE HSIQYSGEVR  RFNTANDDNV TQVRAFVNV LNEEQKRLC ENIAGHLKDA QIFIQKKAVK  NFTEVHPDYG SHIQALLDKY NAEKPKNAIH TQVQSGSHLA AREKANL</p>
<b>Source</b>	Baculovirus
<b>Target Names</b>	CAT
<b>Protein Names</b>	Recommended name: Catalase EC= 1.11.1.6
<b>Expression Region</b>	2-527
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	<p>This gene encodes catalase, a key antioxidant enzyme in the bodies defense against oxidative stress. Catalase is a heme enzyme that is present in the peroxisome of nearly all aerobic cells. Catalase converts the reactive oxygen species hydrogen peroxide to water and oxygen and thereby mitigates the toxic effects of hydrogen peroxide. Oxidative stress is hypothesized to play a role in the development of many chronic or late-onset diseases such as diabetes, asthma, Alzheimer s disease, systemic lupus erythematosus, rheumatoid arthritis, and cancers. Polymorphisms in this gene have been associated with decreases in catalase activity but, to date, acatalasemia is the only disease known to be caused by this gene.</p>
<b>Reconstitution</b>	<p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final</p>



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.