



Recombinant Type-2 restriction enzyme EcoRII (ecoRIIR)

Product Code	CSB-MP318560ENL
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P14633
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	MLMSVFHNWL LEIACENYFV YIKRLSANDT GATGGHQVGL YIPSGIVEKL FPSINHTREL NPSVFLTAHV SSHDCPDSEA RAIYNSRHF GKTRNEKRIT RWVEAAHFRI LKITGALTLL AFKLDEQGGD CKEVNIWVCA STDEEDVIET AIGEVI PGAL ISGPAGQILG GLSLQQAPVN HKYILPEDWH LRFPSGSEII QYAASHYVKN SLDPDEQLLD RRRVEYDIFL LVEELHVLDI IRKGFSGSVDE FIALANSVSN RRKS RAGKSL ELHLEHLFIE HGLRHFATQA ITEGNKKPDF LFPSAGAYHD TEFPVENLRM LAVKTTCKDR WRQILNEADK IHQVHLFTLQ EGVSLAQYRE MRESGVRLVV PSSLHKKYPE AVRAELMTLG AFIAELTGLY ADIP
Source	Mammalian cell
Target Names	ecoRIIR
Protein Names	Recommended name: Type-2 restriction enzyme EcoRII Short name= R.EcoRII EC= 3.1.21.4 Alternative name(s): Endonuclease EcoRII Type II restriction enzyme EcoRII
Expression Region	1-404
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.