



# Recombinant Type-2 restriction enzyme EcoRII (ecoRIIR)

<b>Product Code</b>	CSB-MP318560ENL
<b>Storage</b>	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
<b>Uniprot No.</b>	P14633
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MLMSVFHNWL LEIACENYFV YIKRLSANDT GATGGHQVGL YIPSGIVEKL FPSINHTREL NPSVFLTAHV SSHDCPDSEA RAIYYNSRHF GKTRNEKRIT RWVEAAHFRI LKITGALTLL AFKLDEQGGD CKEVNIWVCA STDEEDVIET AIGEVIPGAL ISGPAGQILG GLSLQQAPVN HKYILPEDWH LRFPSGSEII QYAASHYVKN SLDPDEQLLD RRRVEYDIFL LVEELHVLDI IRKGFSGSVDE FIALANSVSN RRKS RAGKSL ELHLEHLFIE HGLRHFATQA ITEGNKKPDF LFPSAGAYHD TEFVENLRM LAVKTTCKDR WRQILNEADK IHQVHLFTLQ EGVSLAQYRE MRESGVRLVV PSSLHKKYPE AVRAELMTLG AFIAELTGLY ADIP
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ecoRIIR
<b>Protein Names</b>	Recommended name: Type-2 restriction enzyme EcoRII Short name= R.EcoRII EC= 3.1.21.4 Alternative name(s): Endonuclease EcoRII Type II restriction enzyme EcoRII
<b>Expression Region</b>	1-404
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.