



Recombinant Rhodobacter sphaeroides Modification methylase Rsrl (rsrIM)

Product Code	CSB-MP318583RLE
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P14751
Product Type	Recombinant Protein
Immunogen Species	Rhodobacter sphaeroides (Rhodopseudomonas sphaeroides)
Purity	≥85% (SDS-PAGE)
Sequence	MANRSHHNAG HRAMNALRKS GQKHSSSQL GSSEIGTTRH VYDVCDCLDT LAKLPDDSVQ LIICDPYNI MLADWDDHMD YIGWAKRWLA EAERVLSPTG SIAIFGGLQY QGEAGSGDLI SIISHMRQNS KMLLANLIW NYPNGMSAQR FFANRHEEIA WFAKTKKYFF DLDVAVREPYD EETKAAYMKD KRLNPESVEK GRNPTNVWRM SRLNGNSLER VGHPTQKPAA VIERLVRALS HPGSTVLDFD AGSGVTARVA IQEGRNSICT DAAPVFKEYY QKQLTFLQDD GLIDKARSYE IVEGAANFGA ALQRGDVAS
Source	Mammalian cell
Target Names	rsrIM
Protein Names	Recommended name: Modification methylase Rsrl Short name= M.Rsrl EC= 2.1.1.72 Alternative name(s): Adenine-specific methyltransferase Rsrl
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.