



# Recombinant *Oryza sativa* subsp. japonica Barley B recombinant-like protein B (Os10g0115200, LOC\_Os10g02584)

<b>Product Code</b>	CSB-MP318680OFG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 mo
<b>Uniprot No.</b>	P0DH89
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDDDASMSIR WGGFFESPAR NLGLQLMSSV PAERDTKQLL SGSPFLHHQH QQHVPHHHHQ PHHPRDCGAN GNANGGAMPP PPATEAPPSM PMNFARSDMW MHPQQQQQHH HPREHKALHN LTVGHGSSHI AHHDPVGYGM IPGHTLQMM QQQTEPQLQP PPPPQQPKEE CISSPLIEEN VPVIDEPPPP KKRQQGRQPK VPRAKKPKKS AAPREDGAPP NAPAPRRRGP RKNIGMNING IDLDLSRIPT PICSCTGAPQ QCYRWGAGGW QSACCTTTIS TYLPMSTKR RGARIAGRKM SHGAFKKVLE KLAGEGYNLN NPIDLKTFWA KHGTNKFVTI R
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Os10g0115200
<b>Protein Names</b>	Recommended name: Barley B recombinant-like protein B Short name= BBR-like protein B Alternative name(s): GAGA-binding transcriptional activator BBR-B
<b>Expression Region</b>	1-341
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.