



Recombinant Arabidopsis thaliana F-box protein At5g36820 (At5g36820)

Product Code	CSB-EP318681DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 mo
Uniprot No.	P0DI01
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAMSNLPRDL LEEVLSRVPV KSIAAVRSTC KNWNSLTYGQ SFTKKLYGKT MATKEKEFLV VMTMDLEVYL MRVNLHGIHK DDNNVKSSIM QKAKLIRLND DRVRVDDICK VFHCDGLLLC ITIGIRLVVC NPYCGQTRCI KTRRDYHITD NYALGHEKMK NSPLRNYKIL VFHDKSFLQN SWFEIYNFNS DSWKVLVYFTC DWKLPFSQLV VSLKGNTYWF AREMYIHGPR IDLPDFLICF DFTTERFGPR LHLPFHSRCV DTVTLASVRE EQLAVLFQDS KTLILEVWIT TKIEPNAVSW SSKVFLEVNM SPLTGFQFNR SFGSFFIVEE KNVVVVIPIKG GHFKRNLAYI IGKDEYFKEV DLGVPSSYIY FSPHVCSYVP SLVQIKKDAQ VMLQHHNVSA EKHHEFCASL
Source	E.coli
Target Names	At5g36820
Protein Names	Recommended name: F-box protein At5g36820
Expression Region	1-410
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.