



Recombinant Drosophila melanogaster Protein twist (twi)

Product Code	CSB-EP318774DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P10627
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MMSARSVSPK VLLDISYKPT LPNIMELQNN VIKLIQVEQQ AYMQSGYQLQ HQQQHLHSHQ HHQQHHQQQH AQYAPLPSEY AAYGITELED TDYNIPSNEV LSTSSNQSAQ STSLELNNNN TSSNTNSSGN NPSGFDGQAS SGSSWNEHGK RARSSGDYDC QTGGSLVMQP EHKKLIHQQQ QQQQQHQQQI YVDYLPTTV D EVASAQSCPG VQSTCTSPQS HFDFFDEELP EHKAQVFLPL YNNQQQQSQQ LQQQQPHQQS HAQMHFQNAV RQSFEGYEPA NSLNGSAYSS SDRDDMEYAR HNALSSVSDL NGGVMSPACL ADDGSAGSLL DGSDAGGKAF RKPRRRLKRK PSKTEETDEF SNQRVMANVR ERQRTQSLND AFKSLQQIIP TLP SDKLSKI QTLKLATRYI DFLCRMLSSS DISLLKALEA QGSPSAYGSA SLLSAAANG AEADLKCLRK ANGAPIIPPE KLSYLFGVWR MEGDAQHQKA
Source	E.coli
Target Names	twi
Protein Names	Recommended name: Protein twist
Expression Region	1-490
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.