



Recombinant Human Semenogelin-1 (SEMG1)

Product Code	CSB-YP021001HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P04279
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>QKGGSKGRLPSEFSQFPHGQKQGQHYSGQKQKQQTESKGSFSIQYTYHVDAN DHDQSRKSQQYDLNALHKTTKSQRHLGGSQQLLNKQEGRDHDKSKGHFHR VVIHHKGGKAHRGTQNPSSQDQGNPSGKGISSQYSNTEERLWVHGLSKEQT SVSGAQKGRKQGGSSSYVLQTEELVANKQQRETKNSHQKNGHYQNVVEV REEHSSKVQTSLCPAHQDKLQHGSKDIFSTQDELLVYNKNQHQTKNLNQDQQ HGRKANKISYQSSSTEERLHYGENGVQKDVSRSIYSQTEKLVAGKSQIQAP NPKQEPWHGENAKGESGQSTNREQDLLSHEQKGRHQHGSHGGLDIVIIEQED DSDRHLAQLNNDRNPLFT</p>
Source	Yeast
Target Names	SEMG1
Protein Names	Recommended name: Semenogelin-1 Alternative name(s): Semenogelin I Short name= SGIC Cleaved into the following 3 chains: 1. Alpha-inhibin-92 2. Alpha-inhibin-31 3. Seminal basic protein
Expression Region	24-402
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is the predominant protein in semen. The encoded secreted protein is involved in the formation of a gel matrix that encases ejaculated spermatozoa. The prostate-specific antigen (PSA) protease processes this protein into smaller peptides, with each possibly having a separate function. The proteolysis process breaks down the gel matrix and allows the spermatozoa to move more freely.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.