



# Recombinant Mouse Myelin basic protein (Mbp)

<b>Product Code</b>	CSB-MP013551MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P04370
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGNHSGKREL SAEKASKDGE IHRGEAGKKR SVGKLSQTAS ESDSVFGEAD AIQNNGTSAE DTAVTDSKHT ADPKNNWQGA HPADPGNRPH LIRLFSRDAP GREDNTFKDR PSEDELQTI QEDPTAASGG LDVMASQKRP SQRSKYLATA STMDHARHGF LPRHRDTGIL DSIGRFFSGD RGAPKRGSGK DSHTRTTHYG SLPQKSQHGR TQDENPVVHF FKNIVTPRTP PPSQGGKGRD SRSGSPMARR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Mbp
<b>Protein Names</b>	Recommended name: Myelin basic protein Short name= MBP Alternative name(s): Myelin A1 protein
<b>Expression Region</b>	1-250
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	The protein encoded by the classic MBP gene is a major constituent of the myelin sheath of oligodendrocytes and Schwann cells in the nervous system. However, MBP-related transcripts are also present in the bone marrow and the immune system. These mRNAs arise from the long MBP gene (otherwise called Golli-MBP ) that contains 3 additional exons located upstream of the classic MBP exons. Alternative splicing from the Golli and the MBP transcription start sites gives rise to 2 sets of MBP-related transcripts and gene products. The Golli mRNAs contain 3 exons unique to Golli-MBP, spliced in-frame to 1 or more MBP exons. They encode hybrid proteins that have N-terminal Golli aa sequence linked to MBP aa sequence. The second family of transcripts contain only MBP exons and produce the well characterized myelin basic proteins. This complex gene structure is conserved among species suggesting that the MBP transcription unit is an integral part of the Golli transcription unit and that this arrangement is important for the function and/or regulation of these genes.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.