



Recombinant *Drosophila melanogaster* Protein terminus (term)

Product Code	CSB-BP318861DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P11455
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MFDEPIISSF VFTNEETTQT FHHQWFSQQG LHECATCYSS IDADEPPSQH WLRGGEASQG LHLTKQQQAV MDIIEARQIE TFFFCDESSK DKLEHFMGET CARGIPELLR WMFQNNTVAV EFNLACYVNA MDQVLIFQSG SLRADHHYDV DESVGVVYEM LMQRIENYLN CSSEYGMAEC SITRLKVQVK RIRVEADGQS ADSSVFALPL QLQEEGLTA TTGCSTSEAE LASLRSAYLK HFRECNGYFP PNMRVNLYGL QQCKTTKELY VVPYHISLTL QQLPNKNFLI LNNIMGQFQR LHELSTPVNS IERDQTSSPL KDLHCRRCRT QFSRRSKLHI HQKLRCGQDF SVDSMHADIV EIYEQCLPIS RSVFQHACYG ITKPKTMMRK GQFVPIECDW RSESSVKVQH GPCVVISNAQ HNSPCKFY
Source	Baculovirus
Target Names	term
Protein Names	Recommended name: Protein terminus
Expression Region	1-428
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.