



Recombinant Escherichia coli Uncharacterized protein yhaC (yhaC)

Product Code	CSB-BP318899ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P11864
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MFPVSSIGND ISSDLVRRKM NDLPESPTGN NLEALAPGIE KLKQTSIEMV TLLNLTLPQGG KCIITGDFQK ELAYLQNVIL YNVSSLRLDF LGYNAQIIQR SDNTCELTIN EPLKNQEIST GNININCPK DIYNEIRRLN VIFSCGTGDI VDLSSDLRN VLDYDFTD KHMANTILNP FKLNSTNFTN ANMFQVNFVS STQNATISWD YLLKITPVLISISDMYSEEK IKFVESCLNE PGDITEEQLK IMRFAIKSI PRATLTDKLE NELTKEIYKS SSKIINCLNR IKLTEMKEFS SEKIYDYIDI IIEDYENTKE NAYLVVPQIN YTMDLNIEDS SSEELLSNT LEKDENSEPDN GFEVGEYNTY EAYNSEKQYF TREDYTYDYD LLNAI
Source	Baculovirus
Target Names	yhaC
Protein Names	Recommended name: Uncharacterized protein yhaC Alternative name(s): ORF B' ORFX
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.