



Recombinant Type III restriction-modification system EcoP15I enzyme mod (mod)

Product Code	CSB-MP318953ENL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P12364
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	≥85% (SDS-PAGE)
Sequence	MKKETIFSEV ETANSKQLAV LKANFPQCFD KNGAFIQEKL LEIRASEVE LSKESYSLNW LGKSYARLLA NLPPKTLLE DKTHNQREEN KNSQHLLIKG DNLEVLKHMV NAYA EKVKMI YIDPPYNTGK DGFVYNDDRK FTPEQLSELA GIDLDQAKRI LEFTTKGSSS HSAWLTFIYP RLYIARELMR EDGTIFISID HNEFSQLKLV CDEIFGEQNH VGDLVWKNAT DNNPSNIAVE HEYIIVYTKN KEQLISEWKS NISDVKNLLV NIGEEFASKY TGNELQEKYT QWFREHRSEL WPLDRYKYID KDGIYTGSSQS VHNPGKEGYR YDIIHPKTKK PCKQQPLMGY RFPLDTMDRL LSEEKIFGD DENKIIELKV YAKDYKQKLS SVIHLDRVA TNELKELFPE MTQPFTNAKT IKLVEDLISF ACDGEGIVLD FFAGSGTTHA TVFNLNKNK TSYQFITVQL DEPKKDKSDA MKHGYNTIFD LTKERLIRAS KKNRDQGFKV YQLMADFRAK DESELTLSNH TFFDDVVLTP EQYDTLLTTW CLYDGSLLTT PIEDVDLGGY KAHLCDGRLY LIAPNFTSEA LKALLQKVDS DKDFAPNKVV FYGSNFESAK QMELNEALKS YANKKSIELD LVVRN
Source	Mammalian cell
Target Names	mod
Protein Names	Recommended name: Type III restriction-modification system EcoP15I enzyme mod Short name= M.EcoP15I EC= 2.1.1.72 Alternative name(s): EcoP15I methyltransferase
Expression Region	1-645
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.