



Recombinant Drosophila melanogaster cAMP-dependent protein kinase catalytic subunit (Pka-C1)

Product Code	CSB-BP318954DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P12370
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	GNNATTSNK KVDAAEVKE FLEQAKEEFE DKWRRNPTNT AALDDFERIK TLGTGSFGRV MIVQHKPTKD YYAMKILDKQ KVVKLKQVEH TLNEKRILQA IQFPFLVSLR YHFKDNSNLY MVLEYVPGGE MFSHLRKVGR FSEPHSRFYA AQIVLAFEYL HYLDLIYRDL KPENLLIDSQ GYLKVTDFGF AKRVKGRWT LCGTPEYLAP EIILSKGYNK AVDWWALGVL VYEMAAGYPP FFADQPIQIY EKIVSGKVRP PSHFGSDLKD LLRNLLQVDL TKRYGNLKAG VNDIKNQKWF ASTDWIAIFQ KKIEAPFIPR CKGPGDTSNF DDYEEEEALRI SSTEKCAKEF AEF
Source	Baculovirus
Target Names	Pka-C1
Protein Names	Recommended name: cAMP-dependent protein kinase catalytic subunit Short name= PKA C EC= 2.7.11.11
Expression Region	2-353
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.