



Recombinant Enterobacteria phage P4 Protein beta (Beta)

Product Code	CSB-BP319040EDH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P13057
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage P4 (Bacteriophage P4)
Purity	>85% (SDS-PAGE)
Sequence	MKVYFENYSY YPALRTRSAE MTGLNNLSYE NKKKILPLIS LGKWPRSEEI QVSLDKSLEV MSNLPFILDV TKDNSHHCAS SFELLSPENG FKNWIEFCSR NDNIIPVVQM PDSAKLRDIS IQARVLEELK GSIAFRIRNL NTDINKTLTS LVSMNSPENA IVFIDLGYIR GNVSAITAAA INSINQIRTE IPEAISVLA TSFPSSVTNF CRENGQSGYI DVIERELHQN IGGSDVAIYG DHGSIHSVVY DNIIGRYVPR IDIALNDSWY FERRPGMNKE GFIEAAKSIL AEYPHYQRED SWGAMIRNA AIGDIAGMGS PAKWIAVRVN LHLNKQIELS EALQYGFDDH EEDLI
Source	Baculovirus
Target Names	Beta
Protein Names	Recommended name: Protein beta
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.