



Recombinant Human Tetratricopeptide repeat protein 26 (TTC26)

Product Code	CSB-EP025219HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A0AVF1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MMLSRAPAV GRGVQHTDKR KKKGRKIPKL EELLSKRDFE GAITLLEFKR HVGEEEDTN LWIGYCAFHL GDYKRALEEY ENATKEENCN SEVWVNLA CT YFFLGMYKQA EAAGFKASKS RLQNRLLFHL AHKFNDEKKL MSFHQNLQDV TEDQLSLASI HYMRSHYQEA IDIYKRILLD NREYLALNVY VALCYYYKLDY YDVSQEVAVLAV YLQQIPDSTI ALNLKACNHF RLYNGRAAEA ELKSLMDNAS SSFEFAKELI RHNLVVFRRG EGALQVLPPL VDVIPPEARLN LVIYYLRQDD VQEAYNLIKD LEPTTPQEYI LKGVVNAALG QEMGSRDHMK IAQQFFQLVG GSASECDTIP GRQCMASCFE LLKQFDDVLI YLNSFKSYFY NDDIFNFNYA QAKAATGNTS EGEEAFLLIQ SEKMKNYIY LSWLARCIM NKKPRLAWEL YLMETSGES FSLLQLIAND CYKMGQFYYS AKAFDVLRL DPNPEYWEGK RGACVGIFQM IAGREPKET LREVLHLLRS TGNTQVEYMI RIMKKWAKEN RVSI
Source	E.coli
Target Names	TTC26
Protein Names	Recommended name: Tetratricopeptide repeat protein 26 Short name= TPR repeat protein 26
Expression Region	1-554
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.