



# Recombinant Enterobacteria phage T4 Fibrin (wac)

<b>Product Code</b>	CSB-BP319157EDZ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P10104
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Enterobacteria phage T4 (Bacteriophage T4)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	TDIVLNDLP FVDGPPAEGQ SRISWIKNGE EILGADTQYG SEGSMNRPTV SVLRNVEVLD KNIGILKTSL ETANSDIKI QGILDVSGDI EALAQIGINK KDISDLKTLT SEHTEILNGT NNTVDSILAD IGPFNAEANS VYRTIRNDLL WIKRELGQYT GQDINGLPVV GNPSSGMKHR IINNTDVITS QGIRLSELET KFIESDVGSL TIEVGNLREE LGPKPPSFSQ NVYSRLNEID TKQTTVESDI SAIKTSIGYP GNNSIITSVN TNTDNIASIN LELNQSGGIK QRLTVIETSI GSDDIPSSIK GQIKDNTTSI ESLNGIVGEN TSSGLRANVS WLNQIVGTDS SGGQPSPPGS LLNRVSTIET SVSGLNNAVQ NLQVEIGNNS AGIKGQVVAL NTLVNGTNPV GSTVEERGLT NSIKANETNI ASVTQEVNTA KGNISLQGD VQALQEAGYI PEAPRDGQAY VRKDGWVFL STFLSPA
<b>Source</b>	Baculovirus
<b>Target Names</b>	wac
<b>Protein Names</b>	Recommended name: Fibrin Alternative name(s): Collar protein Whisker antigen control protein
<b>Expression Region</b>	2-487
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.