



Recombinant Human Mediator of RNA polymerase II transcription subunit 19 (MED19)

Product Code	CSB-MP013652HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A0JLT2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MENFTALFGA QADPPPPPTA LGFGPGKPPP PPPPPAGGGP GTAPPPTAAT APPGADKSGA GCGPFYLMRE LPGSTELTGS TNLITHYNLE QAYNKFCGKK VKEKLSNFLP DLPGMIDLPG SHDNSSLRSL IEKPPILSSS FNPITGTMLA GFRLHTGPLP EQCRLMHIQP PKKKNKHKHKH QSRTQDPVPP ETPSDSDHKK KKKKKEEDPD RKRKKKEKKK KKNRHSPDHP GMGSSQASSS SSLR
Source	Mammalian cell
Target Names	MED19
Protein Names	Recommended name: Mediator of RNA polymerase II transcription subunit 19 Alternative name(s): Lung cancer metastasis-related protein 1 Mediator complex subunit 19
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.