



Recombinant *Saccharomyces cerevisiae*

Uncharacterized protein YER152C (YER152C)

Product Code	CSB-YP319184SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P10356
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MKYKEINFFK GHPSSRLLPR EAVIQATAAI LGPETREYDN DPYNRHPLTY GSDEGALWVR EQICTFLNDQ LFKFENGARS RTRADYLNLN SGASYGMLNI LLQTTLPHNG YTRQAFIITP TYFLINNCFT DAGFKGKMTA INEQGHDSID FESLISALEQ HEAEPQPHST TEMIQGPKLT KKVYRYVMYC IPTFANPSGN TYSLETRRRL IDIARKYDML IITDDVYDIL DYTTPSDELP SPPLRMVHID RSTAPSGEDS FGNTVSNATF SKLIAPGLRF GYHESINANL ARQLSKGGAN VSGGTPSQLN SMIVGEMLRG GAAQRCIAHL RSVYSERATV LTSALKKYMP LGTEIMPLKG GYFTWITLPP AYNAMEISTI LAKKFNVILA DGSNFEVIGD EKNWGQSCFR LSISFLEVDD IDRGIELFGA VCKSHAITNN ITM
Source	Yeast
Target Names	YER152C
Protein Names	Recommended name: Uncharacterized protein YER152C
Expression Region	1-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.