



# Recombinant *Bacillus subtilis* Endoglucanase (eglS)

<b>Product Code</b>	CSB-EP319201BRJ-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P10475
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Bacillus subtilis</i> (strain 168)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	A GTKTPVAKNG QLSIKGTQLV NRDGKAVQLK GISSHGLQWY GEYVNKDSLK WLRDDWGITV FRAAMYADG GYIDNPSVKN KVKEAVEAAK ELGIYVIIDW HILNDGNPNQ NKEKAKEFFK EMSSLYGNTN NVIYEIANEP NGDVNWK RDI KPYAEEVISV IRKNPDNII IVGTGTWSQD VNDAADDQLK DANVMYALHF YAGTHGQFLR DKANYALSKG APIFVTEWGT SDASGNGGVF LDQSREWLKY LDSKTISWVN WNLSDKQESS SALKPGASKT GGWRLSDLSA SGTFVRENIL GTKDSTKDIP ETPSKDKPTQ ENGISVQYRA GDGSMNSNQI RPQLQIKNNG NTTVDLKDVT ARYWYKAKNK GQNFDCDYAQ ICGGNVTHKF VTLHKPKQGA DTYLELGFKN GTLAPGASTG NIQLRLHNDD WSNYAQSGDY SFFKSNTFKT TTKITLYDQG KLIWGTEPN
<b>Source</b>	<i>E. coli</i>
<b>Target Names</b>	eglS
<b>Protein Names</b>	Recommended name: Endoglucanase EC= 3.2.1.4 Alternative name(s): Carboxymethyl-cellulase Short name= CMCase Short name= Cellulase Endo-1,4-beta-glucanase
<b>Expression Region</b>	30-499
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.