



Recombinant Escherichia coli N-acetyl-gamma-glutamyl-phosphate reductase (argC)

Product Code	CSB-MP319298ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P11446
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MLNTLIVGAS GYAGAELVTY VNRHPHMNIT ALTVSAQSND AGKLISDLHP QLKGIVDLPL QPMSDISEFS PGVDVVFLAT AHEVSHDLAP QFLEAGCVVF DLSGAFRVND ATFYEKYYGF THQYPELLEQ AAYGLAEWCG NKLKEANLIA VPGCYPTAAQ LALKPLIDAD LLDLNQWPVI NATSGVSGAG RKA AISNSFC EVSLQPYGVF THRHQPEIAT HLGADVIFTP HLGNFPRGIL ETITCRLKSG VTQAQVAQVL QQAYAHKPLV RLYDKGVPAL KNVVGLPFCD IGFVQGEHL IIVATEDNLL KGAAQAVQC ANIRFGYAET QSLI
Source	Mammalian cell
Target Names	argC
Protein Names	Recommended name: N-acetyl-gamma-glutamyl-phosphate reductase Short name= AGPR EC= 1.2.1.38 Alternative name(s): N-acetyl-glutamate semialdehyde dehydrogenase Short name= NAGSA dehydrogenase
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.