



Recombinant *Saccharomyces cerevisiae* Choline-phosphate cytidyltransferase (PCT1)

Product Code	CSB-BP319514SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P13259
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MANPTTGKSS IRAKLSNSSL SNLFKKNKNK RQREETEEQD NEDKDESKNQ DENKDTQLTP RKRRRLTKEF EEKEARYTNE LPKELRKYRP KGFRFNLPPT DRPIRIYADG VFDFLHGLHM KQLEQCKKAF PNVTLIVGVP SDKITHKLKG LTVLTDKQRC ETLTHCRWVD EVVPNAPWCV TPEFLLEHKI DYVAHDDIPY VSADSDDIYK PIKEMGKFLT TQRTNGVSTS DIITKIIRDY DKYLMRNFAR GATRQELNVS WLKKNLEFK KHINEFRSYF KKNQTNLNNNA SRDLYFEVRE ILLKKTGKK LYSKLIGNEL KKQNQRQRKQ NFLDDPFTRK LIREASPATE FANEFTGENS TAKSPDDNGN LFSQEDDED TNSNNTNTNSD SDSNTNSTPP SEDDDDNDRL TLENLTQKKK QSAN
Source	Baculovirus
Target Names	PCT1
Protein Names	Recommended name: Choline-phosphate cytidyltransferase EC= 2.7.7.15 Alternative name(s): CTP:phosphocholine cytidyltransferase Short name= CCT Short name= CT Phosphorylcholine transferase
Expression Region	1-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.