



Recombinant Escherichia coli RNA polymerase sigma factor rpoS (rpoS)

Product Code	CSB-BP319548ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P13445
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSQNTLVKVDH LNEDAEFDEN GVEVFDEKAL VEQEPSDNDL AEEELLSQGA TQRVLDTQL YLGEIGYSPL LTAEDEVYFA RRALRGDVAS RRRMIESNLR LVVKIARRYG NRGLALLDLI EEGNLGLIRA VEKFDPERGF RFSTYATWWI RQTIERAIMN QTRTIRLPIH IVKELNVYLR TARELSHKLD HEPsAEEIAE QLDKPVDDVS RMLRLNERIT SVDTPGGDS EKALLDILAD EKENGPEDTT QDDDMKQSIV KWL FELNAKQ REVLARRFGL LGYEAATLED VGREIGLTRE RVRQIQVEGL RRLREILQTQ GLNIEALFRE
Source	Baculovirus
Target Names	rpoS
Protein Names	Recommended name: RNA polymerase sigma factor rpoS Alternative name(s): Sigma-38
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.