



# Recombinant Pisum sativum Vicilin

<b>Product Code</b>	CSB-MP319617EWE
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P13918
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pisum sativum (Garden pea)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	RS DPQNPFIFKS NKFQTLFENE NGHIRLLQKF DQRSKIFENL QNYRLLLEYKS KPHTIFLPQH TDADYILVVL SGKAILTVLK PDDRNSFNLE RGDTIKLPAG TIAYLVNRDD NEELRVLDLA IPVNRPGQLQ SFLLSGNQNQ QNYLSGFSKN ILEASFNTDY EEIEKVLLEE HEKETQHRRS LKDKRQQSQE ENVIVKLSRG QIEELSKNAK STSKKSVSSE SEPFNLRSG PIYSNEFGKF FEITPEKNPQ LQDLDFVNS VEIKEGSLLL PHYNSRAIVI VTVNEGKGDF ELVGQRNENQ QEQRKEDDEE EEQGEEEINK QVQNYKAKLS SGDVFVIPAG HPVAVKASSN LDLLGFGINA ENNQRNFLAG DEDNVISIQI RPVKELAFPG SAQEVDRILE NQKQSHFADA QPQQRERGSR ETRDRLSSV
<b>Source</b>	Mammalian cell
<b>Protein Names</b>	Recommended name: Vicilin
<b>Expression Region</b>	29-459
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.