



Recombinant Human Fez family zinc finger protein 1 (FEZF1)

Product Code	CSB-MP008602HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A0PJY2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDSSCHNATT KMLATAPARG NMMSTSKPLA FSIERIMART PEPKALPVP FLQGALPKGE PKHSLHLNSS IPCMIPFVPV AYDTSPKAGV TGSEPRKASL EAPAAPAAVP SAPAFSCSDL LNCALSLKGD LARDALPLQQ YKLVRPRVVN HSSFHAMGAL CYLNRGDGPC HPAAGVNIHP VASYFLSSPL HPQPKTYLAE RNKLVVPAVE KYPGVAFKD LSQAQLQHYM KESAQLLSEK IAFKTSDFSR GSPNAKPKVF TCEVCGKVFN AHYNLTRHMP VHTGARPFVC KVCCKGFRQA STLCRHKIIH TQEKPHKCNQ CGKAFNRSST LNTHTRIHAG YKPFVCEFCG KGFHQKGNKY NHKLTHSGEK QFKCNICNKA FHQVYNLTFH MHTHNDKKPF TCPTCGKGFC RNFDLKKHVR KLHDSSLGLA RTPAGEPGTE PPPPLPQQPP MTLPPLQPPL PTPGPLQPGL HQGHQ
Source	Mammalian cell
Target Names	FEZF1
Protein Names	Recommended name: Fez family zinc finger protein 1 Alternative name(s): Zinc finger protein 312B
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.