



# Recombinant Human Fez family zinc finger protein 1 (FEZF1)

<b>Product Code</b>	CSB-YP008602HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A0PJY2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDSSCHNATT KMLATAPARG NMMSTSKPLA FSIERIMART PEPKALPVP FLQGALPKGE PKHSLHLNSS IPCMIPFVPV AYDTSPKAGV TGSEPRKASL EAPAAPAAVP SAPAFSCSDL LNCALSLKGD LARDALPLQQ YKLVRPRVVN HSSFHAMGAL CYLNRGDGPC HPAAGVNIHP VASYFLSSPL HPQPKTYLAE RNKLVVPAVE KYPGVAFKD LSQAQLQHVM KESAQLLSEK IAFKTSDFSR GSPNAKPKVF TCEVCGKVFN AHYNLTRHMP VHTGARPFVC KVCCKGFRQA STLCRHKIIH TQEKPHKCNQ CGKAFNRSST LNTHTRIHAG YKPFVCEFCG KGFHQKGNKY NHKLTHSGEK QFKCNICNKA FHQVYNLTFH MHTHNDKKPF TCPTCGKGFC RNFDLKKHVR KLHDSSLGLA RTPAGEPGTE PPPPLPQQPP MTLPPLQPPL PTPGPLQPGL HQGHQ
<b>Source</b>	Yeast
<b>Target Names</b>	FEZF1
<b>Protein Names</b>	Recommended name: Fez family zinc finger protein 1 Alternative name(s): Zinc finger protein 312B
<b>Expression Region</b>	1-475
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.