



Recombinant *Drosophila melanogaster* Accessory gland-specific peptide 26Aa (Acp26Aa)

Product Code	CSB-EP319885DLU-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P10333
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	SC DSEQQLDSAM HLKSDSTKSA SLKNVAPKND ETQAKIAKDD VALKDAKKGD YIMDIDISDL PLDDYPINRS KSLKSSSIDL NNIPFNKGLD DFPAKEKNQG SNQSALKALQ QRLLETEQNNS LLLRNHSIYL MKEIEARKTD IIKVRQLNLD LELELNTVNR RLLELNGQLQ NTRKSTKPCK KRSSKDSAPP AANQFQEANV RNTYRNKYL LKELSQKIN NEIAKVATDV PTETNPSQGN LPTL
Source	<i>E.coli</i>
Target Names	Acp26Aa
Protein Names	Recommended name: Accessory gland-specific peptide 26Aa Alternative name(s): Male accessory gland secretory protein 355A
Expression Region	19-264
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.