



Recombinant rRNA adenine N-6-methyltransferase (ermBC)

Product Code	CSB-MP319929ENL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P10738
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	MNKNIKYSQN FLTSEKVLNQ IIKQLNLKET DTVYEIGTGK GHLTTKLAKI SKQVTSIELD SHLFNLSSEK LKSNTRVTLI HQDILQFQFP NKQRYKIVGN IPYHLSTQII KKVVFESHAS DIYLIVEEGF YKRTLDIHRT LGLLLHTQVS IQQLLKLP AE CFHPKPRVNS VLIKLRHTT DVPDKYWKLY TYFVSKWVNR EYRQLFTKNQ FHQAMKHAKV NNLSTVTYEQ VLSIFNSYLL FNGRK
Source	Mammalian cell
Target Names	ermBC
Protein Names	Recommended name: rRNA adenine N-6-methyltransferase EC= 2.1.1.184 Alternative name(s): Erythromycin resistance protein Macrolide-lincosamide-streptogramin B resistance protein
Expression Region	1-245
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.