



# Recombinant Rat Coiled-coil domain-containing protein 61 (Ccdc61)

<b>Product Code</b>	CSB-EP004722RA-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A0JPP8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEQPAGLQVD YMFRGVEHAV RVVVSQGQVLE LEVEDRMTAD QWRGEFDANF IEDLTHKTGN FKQFSIFCNM LESALTQSSE SVTLDLLTYT DLESLRSRKL GGRPGPCPRS AQLNSKRYLI LIYSVEFDRI HYPLPLPYQG KPDPVVLLQGI IRSLKEELGH LRGMNGGQDA RETEIWHLRE QVTRLASEKR ELEAQLGRSR EEALAGRAAR QEAESLRGLV RGLELELRQE RGLGGRAAGR RSQDCRRLAK ELEEKASER NLRARLKTIN CELAMYRRGR RTLPAGARED RALSSRERST SRGRTATRSS SRESNRGARS HGRPAHPSPS PTGSRVPRFD PTAFVKAKEK KQREIRMKRQ QQQQQQRNRM GSGGSGDGPS VSWSHQTRPP AAVTGRGDAA NRSRNRSSSV DSFRSRCSSV SSCSELEDFS QSVSKSRRCR GRGKPPSPIP WSGSKTKSTT RERNNHQRHL ASSGAWVPIK EYSSDYQGAD MAEIDARLKA LQEYMNRLDT RS
<b>Source</b>	E.coli
<b>Target Names</b>	Ccdc61
<b>Protein Names</b>	Recommended name: Coiled-coil domain-containing protein 61
<b>Expression Region</b>	1-512
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.