



Recombinant Transposase for insertion sequence element IS421

Product Code	CSB-MP320045ENL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P11901
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	≥85% (SDS-PAGE)
Sequence	MNYSHDNWSA ILAHIGKPEE LDTSARNAGA LTRRREIRDA ATLLRLGLAY GPGGMSSLRE VTAWAQLHDV ATLSDVALLK RLRNAADWFG ILAAQTLAVR AAVTGCTSGK RLRLVDGTAI SAPGGGSAEW RLHMGYDPHT CQFTDFELTD SRDAERLDRF AQTADEIRIA DRGFGSRPEC IRSLAFGEAD YIVRVHWRGL RWLTAEGMRF DMMGFLRGLD CGKNGETTVM IGNSGNKKAG APFPARLIAV SLPPEKALIS KTRLLSENRR KGRVVQAETL EAAGHVLLLT SLPEDEYSAE QVADCYRLRW QIELAFKRLK SLLHLDALRA KEPELAKAWI FANLLAAFLI DDIIQPSLDF PPRSAGSEKK N
Source	Mammalian cell
Protein Names	Recommended name: Transposase for insertion sequence element IS421
Expression Region	1-371
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.