



Recombinant *Oryza sativa* subsp. *japonica* 30S ribosomal protein S12, chloroplastic (rps12-A)

Product Code	CSB-EP320067OFG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P12149
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MPTVKQLIRN ARQPIRNARK SAALKGCPQR RGTCARVYTI NPKKPNSALR KVARVRLTSG FEITAYIPGI GHNLQEHSVV LVRGGRVKDL PGVRYRIIRG ALDAVAVKNR QQGRSKYGVK KPKK
Source	E.coli
Target Names	rps12-A;
Protein Names	Recommended name: 30S ribosomal protein S12, chloroplastic
Expression Region	1-124
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.