



Recombinant *Marchantia polymorpha* Maturase K (matK)

Product Code	CSB-MP320069MCB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P12174
Product Type	Recombinant Protein
Immunogen Species	<i>Marchantia polymorpha</i> (Liverwort) (<i>Marchantia aquatica</i>)
Purity	>85% (SDS-PAGE)
Sequence	MEHRIYNSNY FLDITIPYFF HPEILIRIFR RHIQDIPFLH FLRTLTYKNK CLNILNIENF FYLKKNQFFC FLWNFYIYEF EYLLNDIWEK FYKFESVFFW NFIDKTNSIK KIKHILKSK KPIEKKIVKK ISSIHYIRYK NNLITLNDR NILILENWKD FFLIFWQKYF NVWFKSSRIL IQNFYKNSFS FLGYMFRIES QIILIQIQL NLLRNVNLK KEFCSIIPVI PLIRLLAKEK FCDVLGRPLC KLSWTTLSDN EIFERFDQII KHIFSYYSGC INKKGLYQLQ YIFRFSCAKT LACKHKSTIR TVWKKYGSNL LTSSIFFNKT KLISLNFNSK NPYKKNFWYL NIIQVNYLAH SLQSKLLKE
Source	Mammalian cell
Target Names	matK
Protein Names	Recommended name: Maturase K Alternative name(s): Intron maturase
Expression Region	1-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.