



Recombinant Human immunodeficiency virus type 1 group M subtype B Truncated surface protein (env)

Product Code	CSB-MP320112HKY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P12490
Product Type	Recombinant Protein
Immunogen Species	Human immunodeficiency virus type 1 group M subtype B (isolate NY5) (HIV-1)
Purity	≥85% (SDS-PAGE)
Sequence	QLWVTVYYG VPVWKEATTT LFCASDAKAY DTEVHNWVAT HACVPTDPNP QEVVLQNVTE NFNMWKNNTV EQMHEDIISL WDQSLKPCVK STPLCVTLNC TDLTNATYAN GSSEERGEIR NCSFNVTTII RNKIQKEYAL FYRLDIVPID KDNTSYTLIN CNTSVITQAC PKVSFEPIPI HYCAPAGFAI LKCNDKKFNG TGPCTNVSTV QCTHGIKPVV STQLLLNGSL AEGEVVIRSE NFTNNAKTII VQLNKSVEIN CTRPNNNTKK GIAIGPGRTL YAREKIIGDI RQAHCNISKA KWNDTLKQIV TKLKEQFRNK TIVFNQSSGG DPEIVMHSFN CGGEFFYCKT TQLFNSTWLF NSTWNDTERS DNNETIIIPC RIKQIINSGR K
Source	Mammalian cell
Target Names	env
Protein Names	Recommended name: Truncated surface protein Short name= SU Alternative name(s): Glycoprotein 120 Short name= gp120
Expression Region	32-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.