



# Recombinant Human Shootin-1 (KIAA1598)

<b>Product Code</b>	CSB-EP012273HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A0MZ66
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MNSSDEEKQL QLITSLKEQA IGEYEDLRAE NQKTKEKCDK IRQERDEAVK KLEEFQKISH MVIEEVNFMQ NHLEIEKTCR ESAEALATKL NKENKTLKRI SMLYMAKLGPDVITEEINID DEDSTTDTDG AAETCVSVQC QKQIKELRDQ IVSVQEEKKI LAIELENLKS KLVEVIEEVN KVKQEKTVLN SEVLEQRKVL EKCNRVSMMLA VEEYEEMQVN LELEKDLRKK AESFAQEMFI EQNKLKRQSH LLLQSSIPDQ QLLKALDENA KLTQQLEEER IQHQQKVKEE EQLENETLH KEIHNLKQQL ELLEEDKKEL ELKYQNSEEK ARNLKHSVDE LQKRVNQRSEN SVPPPPPPPP PLPPPPPNPI RSLMSMIRKR SHPSGSGAKK EKATQPETTE EVTDLKRQAV EEMMDRIKKG VHLRPVNQTA RPKTKPESSK GCESAVDELK GILGTLNKST SSRSLKSLDP ENSETELERI LRRRKVTAEA DSSSPTGILA TSEKSMPLV GSVSSVTKTA LNKKTLEAEF NSPSPPTPEP GEGPRKLEGC TSSKVTFQPP SSIGCRKKYI DGEKQAEPVV VLDPVSTHEP QTKDQVAEKD PTQHKEDEGE IQPENKEDSI ENVRETDSSN C
<b>Source</b>	E.coli
<b>Target Names</b>	SHTN1
<b>Protein Names</b>	Recommended name: Shootin-1
<b>Expression Region</b>	1-631
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.