



Recombinant Drosophila melanogaster Vitelline membrane protein Vm26Ab (Vm26Ab)

Product Code	CSB-YP320194DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P13238
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	IQLQPTQ GILIPAPLAE NIRVSRAAYG GYGAAPAAPS YSAPAAPAAQ AYSAPAAPAY SAPAAPAYSA PAAPAYSAPA APAYSAPAAP AYSAPASIPS PPCPKNYLFS CQPSLQPVPC SAPAQSYGSA GAYSQYVPQY AVPFVREL
Source	Yeast
Target Names	Vm26Ab
Protein Names	Recommended name: Vitelline membrane protein Vm26Ab Alternative name(s): Protein SV23 Protein TU-4
Expression Region	24-168
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.