



Recombinant Human Calcium and integrin-binding family member 4 (CIB4)

Product Code	CSB-EP005429HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A0PJX0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGQCLRYQMH WEDLEEYQAL TFLTRNEILC IHDTFLKLCP PGKYYKEATL TMDQVSSLPA LRVNPFDRDI CRVFSHKGMF SFEDVLGMAS VFSEQACPSL KIEYAFRIYD FNENGFIDEE DLQRIILRLI NSDDMSEDLL MDLTNHVLSE SDLDNDNMLS FSEFEHAMAK SPDFMNSFRI HFWGC
Source	E.coli
Target Names	CIB4
Protein Names	Recommended name: Calcium and integrin-binding family member 4
Expression Region	1-185
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.