



Recombinant Escherichia coli Mannosyl-D-glycerate transport/metabolism system repressor mngR (mngR)

Product Code	CSB-YP320258ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P13669
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MGHKPLYRQI ADRIREQIAR GELKPGDALP TESALQTEFG VSRVTVRQAL RQLVEQQILE SIQSGSITYVK EERVNYDIFQ LTSFDEKLSD RHVDTHSEVL IFEVIPADDF LQQQLQITPQ DRVWHVVRVYRQKQKPMALP ETWMPALALFP DLTWQVMENS KYHFIEEVKK MVIDRSEQEI IPLMPTEEMS RLLNISQTKP ILEKVSRYGL VDGRVFEYSR NAFNTDDYKF TLIAQRKSSR
Source	Yeast
Target Names	mngR
Protein Names	Recommended name: Mannosyl-D-glycerate transport/metabolism system repressor mngR Alternative name(s): Fatty acyl-responsive regulator Protein P30
Expression Region	1-240
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.