



Recombinant *Schizosaccharomyces pombe* Serine/threonine-protein phosphatase PP1-1 (dis2)

Product Code	CSB-BP320260SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P13681
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSNPDVDLDS IIDRLLEVRG SRPGRQVQLS EDEIRFLCNK AREIFISQPI LLELEAPLKI CGDIHGQYYD LLRLFYGGF PPEANYLFLG DYVDRGKQSL EVICLLLAYK IKYPENFFIL RGNHECASIN RIYGFYDECK RRYNIKLWKT FTDCFNCLPI AAIIDEKIFT MHGGLSPDLN SMDQIQRIMR PTDVPDTGLL CDLLWSDPDK DLTGWGDNDR GVSFTFGPDV VSRFLHKHDM DLVCRAHQVV EDGYEFFSKR QLVTLFSAPN YCGEFDNAGA MMSVDESLLC SFQILKPAEK KQRYGYQGSS QNWHMTPPRK NKTGNSK
Source	Baculovirus
Target Names	dis2
Protein Names	Recommended name: Serine/threonine-protein phosphatase PP1-1 EC=3.1.3.16
Expression Region	1-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.