



Recombinant Bovine Thioesterase superfamily member 4 (THEM4)

Product Code	CSB-YP023494BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A1A4L1
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	STEK VIHKDWALPN PSWSKDLKLL FDQFMKKCED GSWERLPSYK RRSTQESEDF KTYFLDPKLV EEERLSQAQL FTRGFEDGLG FEYVIFKNND EKRTVCLFQG GPYLQGV PGL LHGGAMATMI DIALGSCTGG AVMTANLNIN FKRPVPLCSV VVINSQLDKL EGRKLF LSCN VRSVDEK TLY SEATGLFIKL DPEKSST
Source	Yeast
Target Names	THEM4
Protein Names	Recommended name: Thioesterase superfamily member 4
Expression Region	37-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Protein kinase B (PKB) is a major downstream target of receptor tyrosine kinases that signal via phosphatidylinositol 3-kinase. Upon cell stimulation, PKB is translocated to the plasma membrane, where it is phosphorylated in the C-terminal regulatory domain. This protein negatively regulates PKB activity by inhibiting phosphorylation. Transcription of this gene is commonly downregulated in glioblastomas.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.