



Recombinant Bovine Protein C-ets-2 (ETS2)

Product Code	CSB-BP007853BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A1A4L6
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	MNDFGIKNMD QVAPVASSYR GTLKRQAAFD TFDGSLLAVF PSLNEEQTLQ EVPTGLDSIS HDSANCELPL LTPCSKAVMS QALKATFSGF KKEQRRLGIP KNPWLWTEQQ VCQWLLWATN EFSLVDVNLQ RFGMTGQVLC NLGKERFLEL APDFVGDILW EHLEQMIKEN QEKNEDQYEE NSHLNSVPHW INSNSLGFGV EQAPYGMQTQ SYPKGGLLDG LCPASSAPST LGPEQDFQMF PKARLNTVSV NYCSVGQDFP AGSLNLLSSA SGKPRDHDSA ETGGDSFESS ESLLQSWNSQ SLLLDVQRVP SFESFEDDCS QSLGLSKPTM SFKDYIQDRS DPVEQGKPMV PAAVLGFTG SGPIQLWQFL LELLSKSCQ SFISWTGDGW EFKLADPDEV ARRWGKRKNK PKMNYEKLSR GLRYYYDKNI IHKTSGKRYV YRFVCDLQNL LGFTPEELHA ILGVQPDTED
Source	Baculovirus
Target Names	ETS2
Protein Names	Recommended name: Protein C-ets-2
Expression Region	1-470
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.