



# Recombinant *Oryza sativa* subsp. *japonica* Glutamine synthetase, chloroplastic (GLN2)

<b>Product Code</b>	CSB-YP320404OFG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P14655
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	VS TETTGVVTRM EQLLNMDTTP FTDKIIAEYI WVGGTGIDLR SKSRTISKPV EDPSELPKWN YDGSSTGQAP GEDSEVILYP QAIFKDPFRG GNNILVMCDT YTPAGEPIPT NKRNRAAQVF SDPKVVSQVP WFGIEQEYTL LQRDVNWPLG WPVGGYPGPQ GPYYCAVGSD KSFGRDISDA HYKACLYAGI NISGTNGEVM PGQWEYQVGP SVGIEAGDHI WISRYILERI TEQAGVVLTL DPKPIQGDWN GAGCHTNYST KSMREDGGFE VIKKAILNLS LRHDLHISAY GEGNERRLTG LHETASIDNF SWGVANRGCS IRVGRDTEAK GKGYLEDRRP ASNMDPYVVT ALLAETTILW EPTLEAEVLA AKKLALKV
<b>Source</b>	Yeast
<b>Target Names</b>	GLN2
<b>Protein Names</b>	Recommended name: Glutamine synthetase, chloroplastic EC= 6.3.1.2 Alternative name(s): Glutamate--ammonia ligase OsGS2 Short name= GS2
<b>Expression Region</b>	49-428
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.