



Recombinant Enterobacteria phage phi80 Repressor protein CI (CI)

Product Code	CSB-EP320431EDL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P14819
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage phi80 (Bacteriophage phi-80)
Purity	≥85% (SDS-PAGE)
Sequence	SSISERIKF LLAREGLKQR DLAEALSTSP QTVNNWIKRD ALSREAAQQI SEKFGYSLDW LLNGEGSPKK DLESNIPPES EWGTVDWDK NTPLPDDEVE VPFLKDIEFA CGDGRVHDED HNGFKLRFSK ATLRRVGANS DGSGVLCFPA SGDSMEPVIP DGATVAVDTG NKRVIDGELY AINQGD LKRI KQLYRKPGGK ILIRSINRDY DDEEADADV EIIGFVFWYS VLR YRR
Source	E.coli
Target Names	CI
Protein Names	Recommended name: Repressor protein CI
Expression Region	2-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.