



# Recombinant *Saccharomyces cerevisiae* Mitochondrial transcription factor 1 (MTF1)

<b>Product Code</b>	CSB-BP320450SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P14908
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SVPIPGIKD ISKLKFFYGF KYLWNPTVYN KIFDKLDLTK TYKHPEELKV LDLYPGVGIQ SAIFYNKYCP RQYSLLEKRS SLYKFLNAKF EGSPQLILKR DPYDWSTYSN LIDEERIFVP EVQSSDHIND KFLTVANVTG EGSEGLIMQW LSCIGNKNWL YRFGKVKMLL WMPSTTARKL LARPGMHSRS KCSVVREAF DTKLIAISDA NELKGFDSQC IEEWDPILFS AAEIWPTKGK PIALVEMDPI DFDFVDNWD YVTRHLMILK RTPLNTVMDS LGHGGQQYFN SRITDKDLLK KCPIDLTNDE FIYLTKLFME WPFKPDILMD FVDMYQTEHS G
<b>Source</b>	Baculovirus
<b>Target Names</b>	MTF1
<b>Protein Names</b>	Recommended name: Mitochondrial transcription factor 1 Alternative name(s): Mitochondrial transcription factor mtTFB Mitochondrial-specificity factor RF1023
<b>Expression Region</b>	2-341
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.