



Recombinant Arabidopsis thaliana F-box/kelch-repeat protein At1g24881 (At1g24881)

Product Code	CSB-MP320519DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0DI03
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MRRRRCDLQP KRTRMCDLQP KRTSMCDLPP KLVGEKILTR IPITSLRAVR STCKLWNALT KDRVLGKAAA QFLGFMTMDS KVCSVRFHLR RSKEEEEDTM DLSIKQVDLL NQVEISRVIYH CDGLLLCVAK DNSRVVWNP YLGQTRWIRP RTESNIGDSY ALGYDINRNH KILRMVQTRN VSVYRYEIID LRSNSWRVLE VTPNGEMDPN HPLYGVSVKG NTYFFAHEDS SSGEIDEDGD IIDLEDFLLC FDFTTETFGL RLPLPFHSTI DATVTLSCVR DQQLAVLYHN EGLHSDDRFT TVEFWTTSI EPNSVSWSKF LTVDMRPLAL TGVRFDNDMG ATFFIDEDEK VAVVFDLDGY LSTESARYHT AFISGKDGFF KPVTLGVAPN VGEP CPRTGH IPTTYRPLV CSSTYLPSLV QVNQQRKRKE RHV
Source	Mammalian cell
Target Names	At1g24881
Protein Names	Recommended name: F-box/kelch-repeat protein At1g24881
Expression Region	1-433
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.