



Recombinant Human Pepsin A-5 (PGA5)

Product Code	CSB-EP320706HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0DJJ9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	VDEQPLEN YLDMEYFGTI GIGTPAQDFT VVFDTGSSNL WVPSVYCSSL ACTNHNRFNP EDSSTYQSTS ETVSITYGTG SMTGILGYDT VQVGGISDTN QIFGLSETEP GSFLYYAPFD GILGLAYPSI SSSGATPVFD NIWNQGLVSO DLFSVYLSAD DKSGSVVIFG GIDSSYYTGS LNWVPVTVEG YWQITVDSIT MNGETIACAE GCQAIVDGTG SLLTGPTSPI ANIQSDIGAS ENSDGMVVS CSAISSLPDI VFTINGVQYP VPPSAYILQS EGSCISGFQG MNVPTESGEL WILGDVFIRQ YFTVFDRANN QVGLAPVA
Source	E.coli
Target Names	PGA5
Protein Names	Recommended name: Pepsin A-5 EC= 3.4.23.1 Alternative name(s): Pepsinogen-5
Expression Region	63-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.