



Recombinant Human Protein FAM170A (FAM170A)

Product Code	CSB-MP008128HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A1A519
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MKRRQKRKHL ENEESQETAE KGGGMSKSQE DALQPGSTRV AKGWSQGVGE VTSTSEYCSC VSSSRKLIHS GIQRIHRDSP QPQSPLAQVQ ERGETPPRSQ HVSLSSYSSY KTCVSSLCVN KEERGMIYY MQVQMNGVA VSWETEETLE SLEKQPRMEE VTLSEVVRVG TPPSDVSTRN LLSDSEPSGE EKEHEERTES DSLPGSPTVE DTPRAKTPDW LVTMENGFRFC MACCRVFTTM EALQEHVQFG IREGFSCHVF HLTMAQLTGN MESESTQDEQ EEENGNEKEE EEKPEAKEEE GQPTTEEDLGL RRSWSQCPGC VFHSPKDRNS
Source	Mammalian cell
Target Names	FAM170A
Protein Names	Recommended name: Protein FAM170A Alternative name(s): Zinc finger domain-containing protein Zinc finger protein ZNFD
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.