



Recombinant Type I restriction enzyme EcoR124II M protein (hsdM)

| | |
|--------------------------|--|
| Product Code | CSB-EP320763ENL |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P10484 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli |
| Purity | >85% (SDS-PAGE) |
| Sequence | MKMTSIQQRA ELHRQIWQIA NDVRGSVDGW DFKQYVLGAL FYRFISENFS SYIEAGDDSI CYAKLDDSVI TDDIKDDAIK TKGYFIYPSQ LFCNVAAKAN TNDRLNADLN SIFVAIESSA YGYNSEADIK GLFADFDTTS NRLGNTVKDK NARLAAVLKG VEGLKLGDFN EHQIDLFGDA YEFLISNYAA NAGKSGGEFF TPQHVSKLIA QLAMHGQTHV NKIYDPAAGS GSLLLQAKKQ FDNHIIIEGF FGQEINHHTY NLARMNMLH NINYDKFDIK LGNTLTEPHF RDEKPFDAIV SNPPYSVKWI GSDDPTLIND ERFAPAGVLA PKSKADFAFV LHALNYLSAK GRAAIVCFPG IFYRGGAEQK IRQYLVNDDY VETVISLAPN LFFGTTIAVN ILVLSKHKTD TNVQFIDASE LFKKETNNNI LTDAHIEQIM QVFASKEDVA HLAKSVAFET VVANDYNLSV SSVYEAKDNR EIIDIAELNA ELKTTVSKID QLRKDIDAIV AEIEGCEVQK |
| Source | E.coli |
| Target Names | hsdM |
| Protein Names | Recommended name: Type I restriction enzyme EcoR124II M protein Short name= M.EcoR124II EC= 2.1.1.72 |
| Expression Region | 1-520 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.